




DentaComfort pillow

The DentaComfort pillow is used by dentists for the positioning of patients in the dental chair.



DentaComfort vacuum pillow adjusts to the patient head shape on the top and the chair on the back and helps patient to relax and lie still in desired position.

The DentaComfort pillow has a range of benefits and features for both dentists and patients.



Easy to give the patient comfortable support in a variety of positions, quickly

Dentist ref:

- My patient was very calm and relaxed throughout the treatment .



Face of the pillow shaped after patient head.



The back of the pillow locking position.



Hyperextension effortless for the patient

Parent:

- Our children feel more secure and enjoy going to the dentist now.



Use your saliva suction for any second stabilizing of the cushion and it is just as quick to reposition.

Patient ref:

- Just fantastic support to my neck, very relaxing and safe.



GERMA

AB GERMA
Box 2127, SE-291 02 KRISTIANSTAD, SWEDEN
Phone: +46 (0)44 12 30 30

Fax: +46-(0)44 10 31 79
E-mail: info@germa.se
Internet: www.germa.se